American Eagle INSTRUMENTS®



Trusted Heritage. Quality Craftsmanship.

YOUNG™

proudly made in the USA

American Eagle INSTRUMENTS®

OUR STORY

From the mountains of Missoula, MT, every American Eagle Instrument begins with a visionary design and US-made steel that is expertly shaped and hand-honed by master craftsmen. Instruments may seem complicated with their many engineered angles, sophisticated metals and ergonomic grips. But like any tool, the magic is in the performance.

Used with pride by generations of dental professionals across the globe, American Eagle Instruments bestow a simple confidence and trusted heritage of US quality craftsmanship.

CUSTOMER SATISFACTION GUARANTEE

If the customer is not completely satisfied within the first 30 days of using their new instrument, we ask they contact their dental dealer to request a credit or replacement.

LIMITED LIFETIME WARRANTY

Except as otherwise noted, American Eagle Instruments offers a lifetime warranty against defects in materials and workmanship and will, at our discretion, replace any product that fails as a result of any such defect.

Each class of instrument has a life expectancy based on its intended use. Due to various contributing factors which could expedite or prolong the effective life of an instrument, it is only reasonable to determine an average life expectancy.

If within the average life expectancy of an instrument our customer experiences corrosion, breakage, or joint failure, American Eagle Instruments will, at our discretion, replace, repair or credit the full value of that instrument to the original owner with proof of purchase.

TABLE OF CONTENTS

4 INSTRUMENT INFORMATION

12 EDUCATION

13 PERIODONTAL

26 QUIK-TIP™ CUSTOMIZATION

35 DIAGNOSTIC

39 RESTORATIVE

48 AUXILIARY PRODUCTS

HOW TO ORDER

Part #

Sec. A	Sec. B	Sec. C	Sec. D
AE	S	H6-7XP	X



When ordering, do not include spaces in the part number.

Section A - Instrument

AE = American Eagle

Section B - Instrument Style

Periodontal/Hygiene

GA = Gracey Access
DG = Double Gracey™
GDP = Gracey Deep Pocket
II = Implant Instrument
G = Gracey
C = Curette
S = Scaler

Diagnostic

EP = Explorer/Probe Combo
EXP = Explorer
M = Mirror
MH = Mirror Handle
P = Probe

Restorative

B = Burnisher
C/P = Condenser/Plugger
CSP = Cement Spatula
EXC = Excavator
HA = Hatchet
MT = Margin Trimmer
CADC = Carver/ Discoid-Cleoid
CAIPC = Carver/Interproximal
C = Chisel
CA = Carver
C = Composite
E = Endodontic
F= File
P = Packing
W = Waxing

Auxiliary

PL = Pliers

Section C - Instrument Classification

QT = Quik-Tip™
TT = Stainless Steel
(Stainless steel instruments require periodic sharpening)
XP = XP® Sharpen-Free Technology
(No sharpening required)

SECTION D - Handle Designation

X = 3/8" EagleLite® Resin Handle
Z = 3/8" EagleLite® Stainless Steel Handle

NOTE: No handle designation for curettes and scalers will indicate an EagleLite® resin handle is requested.

HOW TO ORDER


American Eagle Instruments' products may be purchased through any dealer partner (visit our web-site for a complete list) or by placing your order direct with us and receive your billing from your supplier.


Call: 800-551-5172 or 406-549-7451
Fax: 406-549-7452
Email: customerservice@am-eagle.com

Outside USA, please contact your local American Eagle Instruments dealer.

= American Eagle Instruments Top Seller

INSTRUMENT TIP INFORMATION





CURETTE & SCALER HANDLES

Colored resin handles offer you easier recognition of your specific instrument pattern.

3/8'' - 9.5mm Resin (X)





3/8'' - 9.5mm Stainless Steel (Z)



1-800-551-5172 | Fax: 406-549-7452 | www.am-eagle.com

XP® SHARPEN-FREE TECHNOLOGY

XP SHARPEN-FREE BENEFITS



Saves Time and Money: XP Sharpen-Free instruments last as long as standard instruments, but without the need to sharpen. By eliminating the need for sharpening, time can be used for more productivity.

Maintains Its Original Sharp Edge: These razor sharp edges start and stay sharp longer than any other instrument on the market.

Thinner Blade Design: Since XP Sharpen-Free Instruments are sharpen-free, they are 10% thinner. This allows for easier access of periodontal pockets and interproximal areas.

More Comfort, Less Fatigue: An instrument that is always sharp means only light pressure is needed to "plane" away calculus.

XP SHARPEN-FREE TECHNOLOGY

	Average Lifespan	Maintenance	Ergonomics	Efficacy	Strokes
Conventional Scalers	6-12 months	Sharpen every 1-2 weeks	Abrasive "pop" off calculus	Tip reshaped every sharpening	1,500 
XP Sharpen-Free Scalers	6-12 months	Never sharpen	Gentle shaving stroke	Retains original tip design	15,000 

HOW TO USE XP® SHARPEN-FREE

Holding the Instrument

Hold the instrument using a light, "exploratory scaling" touch.

Shave through calculus using a few light strokes; do not fracture off calculus.

Benefits

- A lighter touch aids in ergonomic hand protection and minimizes the risk of wrist/hand problems.
- A lighter touch increases instrument longevity.
- Less likely to burnish calculus with sharp blade.
- Enhances patient comfort.
- Decreases scaling noise for auditory comfort.
- Less lateral pressure for scaling efficacy.

Guidelines

- XP Sharpen-Free curettes and scalers are designed for scaling and root planing. Avoid using them on margins and overhangs.
- Use less force, a lighter grip and a "shaving" or "exploratory" stroke.
- Do not sharpen your XP Sharpen-Free instruments.
- With proper use and maintenance, you can expect your XP Sharpen-Free instruments to last as long as your stainless steel instruments, without the need to sharpen.



XP SHARPEN-FREE MAINTENANCE

Follow these suggestions to increase the longevity and performance of your instruments.

Care

- Practice care as these instruments have a sharper edge on a thinner blade.
- Avoid fracturing away calculus.
- Avoid use on margins, overhangs, amalgam and composites.

Sterilization

- Use autoclave sterilization and strictly follow unit manufacturer's recommendations.
- Allow for complete drying cycle.
- Do not use dry heat on resin handles.

Storage

- Thoroughly dry before storing.

Cleaning

- Promptly place instruments in ultrasonic cleaner. *Standard recommended cleaning cycle is 20 minutes.*
- Separate XP Sharpen-Free instruments from all other instruments when cleaning/sterilizing. *Combining dissimilar metals is not recommended.*
- Use pH-balanced cleaners.
- Always rinse with demineralized water after ultrasonic cleaning.
- Dry instruments before sterilization.
- Avoid use of chlorine-, chloride-, NaCl- (sodium chloride) and fluoride-based solutions. *Visit our website for a list of compatible solutions.*

STAINLESS STEEL MAINTENANCE

Follow these suggestions to increase the longevity and performance of your instruments.

Care

- Rinse thoroughly with warm water (we recommend using demineralized water, as tap water may contain elements that are detrimental to the instrument).
- Thoroughly dry after cleaning and prior to sterilization.

Sterilization

- Use autoclave sterilization and strictly follow unit manufacturer's recommendations.
- Sterilize hinged instruments in the open position.
- After instrument care and maintenance, thoroughly inspect instruments for any remaining debris or damage prior to use.
- Do not use dry heat on resin handles.

Cleaning

- Stubborn impurities and debris should be manually removed with a soft brush.
- Clean instruments manually, in an ultrasonic cleaning unit, or washer disinfector *(as per unit manufacturer's instructions).*
- Place instruments into cassettes prior to cleaning *(do not overload cassettes).*



IMPLANT INSTRUMENT MAINTENANCE

Follow these suggestions to increase the longevity and performance of your instruments.

Care

- Avoid use on margins, overhangs, amalgam and composites.

Sterilization

- Use autoclave sterilization and strictly follow unit manufacturer's recommendations.
- Allow for complete drying cycle.



Storage

- Thoroughly dry before storing.

PROBE MAINTENANCE

Follow these suggestions to increase the longevity and performance of your instruments.

Cleaning

- Avoid ultrasonic cleaning which may affect the bond of the color markings to the stainless steel.
- Clean the tip with mild liquid soap prior to sterilization.
- Avoid temperatures in excess of 300°F/150°C in the sterilization process, at which point your color probe markings risk becoming discolored.



Sterilization

- Use autoclave sterilization and strictly follow unit manufacturer's recommendations.
- Allow for complete drying cycle.
- Do not use dry heat on resin handles.

1-800-551-5172 | Fax: 406-549-7452 | www.am-eagle.com

INSTRUMENT TESTIMONIALS



" I no longer stress out over my instruments being dull and not having time to sharpen. **American Eagle Instruments have taken that stress off my plate.**
I love having extra time to spend on other things...like my patients."

Angie Stone, RDH, BS
Founder and CEO, HyLife Oral Health Alliance
Author, Amazon Best Selling Book, *Dying From Dirty Teeth*

"American Eagle Instruments, specifically XP Sharpen-Free instruments, have changed how I practice clinically and I will love them forever for that. **My ergonomics are better, my patients are happier and healthier with the lighter shaving stroke,**
and my practice is able to retain patients because of how much more they enjoy their hygiene appointment!"

Andrew Johnston, RDH
Co-host - A Tale of Two Hygienists Podcast





" I decided to give the American Eagle XP Sharpen-Free technology a try and I am so happy that I did! Although there is a learning curve when instrumenting with XP Sharpen-Free, I immediately saw my efficiency improve...
It cut my scaling time nearly in half with my periodontally involved patients and they stay sharper longer than any of the previous instruments I have used in the past."

Jasmin Haley, RDH, MSDH, CDA
Founder of Beyond the Prophy® LLC

CONTINUING EDUCATION

American Eagle Instruments and Young Dental offer accredited continuing education programs through hands on learning and faculty calibration. Visit **www.am-eagle.com** to see a course catalog or to request a live CE course.

HANDS-ON COURSES

Looking to attend an interactive continuing education course? Our clinical education team will bring a hands-on instrumentation course to your study club, component or state meeting for a fun, collaborative learning experience. To schedule a course, contact us at **education@younginnovations.com.**

CLINICAL EDUCATION MANAGER



L. Whitney Howerton
RDH, MDH

Whitney has shown her dedication to the dental hygiene profession through multiple leadership roles on local, state and national levels. She earned a master's degree in dental hygiene from the University of Tennessee Health Science Center. Whitney worked full time as a clinical and didactic instructor teaching a wide variety of courses including local anesthesia, head, neck and dental anatomy, special patient care, and advanced instrumentation. Whitney has been featured in familiar dental hygiene publications such as Access and Dimensions of Dental Hygiene. She is a past recipient of the Esther Wilkins Future Leader Award and a proud member of the ADHA and ADEA. Whitney enjoys speaking professionally in many settings and is eager to share her love for American Eagle Instruments with you!

PERIODONTAL

14 XP® SHARPEN-FREE THIN

16 DOUBLE GRACEY™

18 STANDARD GRACEY

19 GRACEY MESIAL-DISTAL

19 GRACEY ACCESS

20 GRACEY DEEP POCKET

20 CURETTES

22 SCALERS

24 KITS

25 IMPLANT INSTRUMENTS

25 FILES

XP® SHARPEN-FREE THIN

We took the ideal instruments for periodontal maintenance, and made the blades thinner and the handles lighter. This provides you greater tactile sensitivity and increases access into periodontal pockets. Such a thin working end is only possible with XP Sharpen-Free Technology.

AVAILABLE EXCLUSIVELY IN XP SHARPEN-FREE

Dental PRODUCT SHOPPER BEST PRODUCT



M23 XP SHARPEN-FREE THIN SCALER I AESM23TXPX

The M23 XP Sharpen-Free Thin allows for superior supra and interproximal adaptation, and accessibility on all teeth. It also adapts to mesial and distal line angles of all teeth. Its versatility has the potential to replace all posterior and anterior sickle scalers, streamlining preventive and periodontal treatment setups.



BLACKJACK XP SHARPEN-FREE THIN SCALER I AESBJTXPX

Our very popular supragingival access instrument is now thinner, reaching all areas of supragingival calculus like never before. The Blackjack XP Sharpen-Free Thin has the potential to replace all posterior and anterior sickle scalers, including the H6-7, streamlining preventive and periodontal treatment set ups.



BARNHART 5-6 XP SHARPEN-FREE THIN CURETTE I AECB5-6TXPX

The Barnhart XP Sharpen-Free Thin design kicks universal convenience up a notch, making adaptation easier in a wide range of shallow pockets. Curettes are typically too bulky to easily access anterior pockets without tissue damage, particularly at the line angles. The Barnhart XP Sharpen-Free Thin easily adapts in this area, accessing calculus beyond a curette's reach with bulkier designs.

1-800-551-5172 | Fax: 406-549-7452 | www.am-eagle.com



GRACEY ACCESS 1-2 XP SHARPEN-FREE THIN I AEGA1-2TXPX

The Gracey Access 1-2 XP Sharpen-Free Thin curette is useful in pockets greater than 4mm in all anterior areas, as well as buccal and lingual molar furcations.



GRACEY ACCESS 7-8 XP SHARPEN-FREE THIN I AEGA7-8TXPX

The Gracey Access 7-8 XP Sharpen-Free Thin curette is used in pockets greater than 4mm in pre-molar areas, as well as buccal and lingual molar furcations.



GRACEY ACCESS 13-14 XP SHARPEN-FREE THIN I AEGA13-14TXPX

The Gracey Access 13-14 XP Sharpen-Free Thin curette is designed for use in distal surfaces and pockets greater than 4mm in premolar and molar areas. An extended shank improves access to all distal surfaces and pockets of the terminal molars.



GRACEY ACCESS 15-16 XP SHARPEN-FREE THIN I AEGA15-16TXPX


The Gracey Access 15-16 XP Sharpen-Free Thin is a wonderful improvement on the 11-12. The 15-16 shank design improves angulation for better access in all posterior mesial, buccal and lingual areas while easing hand position and reducing operator fatigue. Accessing posterior concavities and furcations is easy and efficient with the 15-16 shank design.

DOUBLE GRACEY™

Truly a universal instrument, Double Gracey instruments have rounded toes and two (2) cutting edges, combining the best features of Gracey and Universal curettes. Advantages of using Double Graceys include easy adaptation to mesial, distal, buccal and lingual surfaces without switching or repositioning the instrument. The mini versions allow easier access into deeper pockets, tighter interproximal areas and furcations. All Double Gracey instruments are sharpen-free with American Eagle Instruments' patented XP® Sharpen-Free Technology.



The face of a Double Gracey is raised, thus creating two slightly sloped faces, approximately 110 degrees, measured from the terminal (lower) shank. The face and blade angles allow the practitioner to use Double Gracey instruments exactly how they would use traditional Gracey curettes — with the terminal shank parallel to the tooth axis.



*All Double Gracey curettes replace Universal curettes.

DOUBLE GRACEY ANTERIOR & MINI ANTERIOR

- Enables clinician to scale from light to heavy calculus.
- Maneuver between buccal and lingual surfaces without switching position.
- XP Sharpen-Free Technology creates less sensitivity on cementum.

DOUBLE GRACEY POSTERIOR & MINI POSTERIOR

- Allows clinician to keep lower shank parallel to tooth surface.
- Maneuver between mesial and distal surfaces without switching position.





DOUBLE GRACEY™ INSTRUMENTS

AVAILABLE EXCLUSIVELY IN XP SHARPEN-FREE



XP SHARPEN-FREE DOUBLE GRACEY ANTERIOR | AEDGAXPX *Replaces Gracey curettes 1-2, 3-4, 5-6, 7-8, and 9-10 along with Universal curettes.*



XP SHARPEN-FREE DOUBLE GRACEY POSTERIOR | AEDGPXPX *Replaces mesial/distal Gracey curettes 11-12, 13-14, and 15-16 along with Universal curettes.*



XP SHARPEN-FREE MINI DOUBLE GRACEY ANTERIOR | AEDGMAXPX *Replaces Access and mini Gracey curettes 1-2, 3-4, 5-6, 7-8, and 9-10 & Universal curettes.*



XP SHARPEN-FREE MINI DOUBLE GRACEY POSTERIOR | AEDGMPXPX *Replaces Access and mini Gracey curettes 11-12, 13-14, and 15-16 & Universal curettes.*

DOUBLE GRACEY KITS

XP SHARPEN-FREE DOUBLE GRACEY KIT | AEDGKXPX

- DOUBLE GRACEY ANTERIOR
- DOUBLE GRACEY POSTERIOR

XP SHARPEN-FREE DOUBLE GRACEY MINI KIT | AEDGMKXPX

- DOUBLE GRACEY MINI ANTERIOR
- DOUBLE GRACEY MINI POSTERIOR

XP SHARPEN-FREE DOUBLE GRACEY 4 KIT* | AEDGK4XPX

- DOUBLE GRACEY ANTERIOR
- DOUBLE GRACEY POSTERIOR
- DOUBLE GRACEY MINI ANTERIOR
- DOUBLE GRACEY MINI POSTERIOR

** Includes stainless steel cassette*

GRACEY HANDLE



3/8'' - 9.5MM EAGLELITE® **COLOR RESIN DOUBLE END (X)**



3/8'' - 9.5MM EAGLELITE **STAINLESS STEEL (Z)**

STANDARD GRACEY

R= RIGID POINT IN **STANDARD GRACEY**



GRACEY 1-2

AEG1-2**X** AEG1-2**Z**

AEG1-2**RZ**



GRACEY 1-2 XP

XP

AEG1-2XP**X**

AEG1-2XP**Z**



GRACEY 3-4

AEG3-4**X** AEG3-4**Z**



GRACEY 3-4 XP

XP

AEG3-4XP**X**

AEG3-4XP**Z**



GRACEY 5-6

AEG5-6**X** AEG5-6**Z**



GRACEY 5-6 XP

XP

AEG5-6XP**X**

AEG5-6XP**Z**



GRACEY 7-8

AEG7-8**X** AEG7-8**Z**

AEG7-8**RZ**



GRACEY 7-8 XP

XP

AEG7-8XP**X**

AEG7-8XP**Z**



GRACEY 9-10

AEG9-10**X** AEG9-10**Z**



GRACEY 9-10 XP

XP

AEG9-10XPX

AEG9-10XPZ



GRACEY 11-12

AEG11-12**X** AEG11-12**Z**

AEG11-12**RZ**



GRACEY 11-12

XP

AEG11-12XP**X**

AEG11-12XP**Z**



GRACEY 13-14

AEG13-14**X** AEG13-14**Z**

AEG13-14**RZ**



GRACEY 13-14 XP

XP


AEG13-14XP**X**

AEG13-14XP**Z**



GRACEY 15-16

AEG15-16**X** AEG15-16**Z**



GRACEY 15-16 XP

XP

AEG15-16XP**X**

AEG15-16XP**Z**



GRACEY 17-18

AEG17-18**X** AEG17-18**Z**




GRACEY 17-18 XP

XP

AEG17-18XP**X**

AEG17-18XP**Z**

STANDARD GRACEY



Gracey curettes are designed to adapt to specific anatomical areas of the teeth for sub-gingival scaling and root planing. The blade angle allows easier insertion into the sulcus while minimizing soft tissue trauma.

GRACEY MESIAL-DISTAL



GRACEY 11-14MD

AEG11-14**X**



GRACEY 11-14 XP

XP

AEG11-14XP**X**

AEG11-14XP**Z**



GRACEY 12-13MD

AEG12-13**X**



GRACEY 12-13 XP

XP

AEG12-13XP**X**

AEG12-13XP**Z**

GRACEY ACCESS *(50% Shorter Blade and Extended Shank)*



GRACEY ACCESS 00-0

AEGA00-0**X**

AEGA00-0**Z**



GRACEY ACCESS 00-0 XP

XP

AEGA00-0XP**X**

AEGA00-0XP**Z**



GRACEY ACCESS 1-2

AEGA1-2**X**



GRACEY ACCESS 1-2 XP

XP THIN

AEGA1-2XP**X**

AEGA1-2XP**Z**

AEGA1-2TXP**X**



GRACEY ACCESS 3-4

AEGA-3-4**X**



GRACEY ACCESS 3-4 XP

XP

AEGA3-4XP**X**

AEGA3-4XP**Z**



GRACEY ACCESS 5-6

AEGA5-6**X**

AEGA5-6**Z**



GRACEY ACCESS 5-6 XP

XP

AEGA5-6XP**X**

AEGA5-6XP**Z**



GRACEY ACCESS 7-8

AEGA7-8**X**

AEGA7-8**Z**




GRACEY ACCESS 7-8 XP

XP THIN

AEGA7-8XP**X**

AEGA7-8TXP**X**



GRACEY ACCESS 11-12

AEGA11-12**X**

AEGA11-12**Z**



GRACEY ACCESS 11-12 XP

XP

AEGA11-12XP**X**

AEGA11-12XP**Z**



GRACEY ACCESS 13-14

AEGA13-14**X**

AEGA13-14**Z**



GRACEY ACCESS 13-14 XP

XP THIN

AEGA13-14XP**X**

AEGA13-14XP**Z**

AEGA13-14TXP**X**



GRACEY ACCESS 15-16

AEGA15-16**X**

AEGA15-16**Z**




GRACEY ACCESS 15-16 XP

XP THIN

AEGA15-16XP**X**

AEGA15-16TXP**X**

GRACEY ACCESS



With the blade half as long and the terminal shank 3mm longer than that of a standard Gracey, the whole blade is easily placed at root surfaces. Narrow pockets and root furcations are also accessible.

GRACEY DEEP POCKET



GRACEY DP 1-2

AEGDP1-2**X**

AEGDP1-2**Z**



GRACEY DP 1-2 XP

XP

AEGDP1-2XP**X**

AEGDP1-2XP**Z**

AVAILABLE EXCLUSIVELY IN XP SHARPEN-FREE



GRACEY DP 3-4 XP

XP

AEGDP3-4XP**X**



GRACEY 5-6

AEGDP5-6**X**



GRACEY 5-6 XP

XP

AEGDP5-6XP**X**

AEGDP5-6XP**Z**



GRACEY DP 7-8

AEGDP7-8**X**

AEGDP7-8**Z**



GRACEY DP 7-8 XP

XP

AEGDP7-8XP**X**



GRACEY DP 11-12

AEGDP11-12**X**

AEGDP11-12**Z**



GRACEY DP 11-12

XP

AEGDP11-12XP**X**

AEGDP11-12XP**Z**



GRACEY DP 13-14

AEGDP13-14**X**

AEGDP13-14**Z**



GRACEY DP 13-14 XP

XP

AEGDP13-14XP**X**


AEGDP13-14XP**Z**



GRACEY DP 15-16

AEGDP15-16**X**

AEGDP15-16**Z**




GRACEY DP 15-16 XP

XP

AEGDP15-16XP**X**

GRACEY DEEP POCKET



With the terminal shank 3mm longer than that of the standard Gracey curette, this instrument can perform scaling and root planing of pockets 5mm or deeper.

CURETTES



BARNHART 1-2

AECB1-2**X**



BARNHART 1-2 XP

XP

AECB1-2XP**X**



BARNHART 5-6

AECB5-6**X**

AECB5-6**Z**



BARNHART 5-6 XP

XP THIN

AECB5-6XP**X**

AECB5-6XP**Z**

AECB5-6XP**Z**

CURETTES *(continued)*



COLUMBIA 13-14

AECC13-14**X**

AECC13-14**Z**




COLUMBIA 13-14 XP

XP

AECC13-14XP**X**

AECC13-14XP**Z**



COLUMBIA 2L-2R

AECC2L-2R**X**



COLUMBIA 2L-2R XP

XP

AECC2L-2RXP**X**

AECC2L-2RXP**Z**



COLUMBIA 4L-4R

AECC4L-4R**X**

AECC4L-4R**Z**



COLUMBIA 4L-4R XP

XP

AECC4L-4R**X**

AECC4L-4R**Z**



LANGER 1-2

AECL1-2**X**

AECL1-2**Z**




LANGER 1-2 XP

XP

AECL1-2XP**X**

AECL1-2XP**Z**



LANGER 3-4

AECL3-4**X**

AECL3-4**Z**



LANGER 3-4

XP

AECL3-4XP**X**

AECL3-4XP**Z**



LANGER 5-6

AECL5-6**X**



LANGER 5-6 XP

XP

AECL5-6XP**X**



MCCALL 13-14S*

AECM13-14S**X**

AECM13-14S**Z**

**scaler-like tip*



MCCALL 13-14S* XP

XP

AECM13-14SXP**X**

**scaler-like tip*



MCCALL 13-14

AECM13-14**X**

AECM13-14**Z**



MCCALL 13-14 XP

XP

AECM13-14XP**X**

AECM13-14XP**Z**



MCCALL 17-18S

AECM17-18S**X**

AECM17-18S**Z**



MCCALL 17-18S XP

XP

AECM17-18SXP**X**

AECM17-18SXP**Z**



MCCALL 17-18

AECM17-18**X**

AECM17-18**Z**



MCCALL 17-18 XP

XP

AECM17-18XP**X**

AECM17-18XP**Z**



RULES 3-4

AECR3-4**X**

AECR3-4**Z**



RULES 3-4 XP

XP

AECR3-4XP**X**

AECR3-4XP**Z**



YOUNGER-GOOD 7-8

AECY7-8**X**

AECY7-8**Z**



YOUNGER-GOOD 7-8 XP

XP

AECY7-8XP**X**

AECY7-8XP**Z**


SCALERS



SCALER 107-108

SCALER 107-108**X**

SCALER 107-108**Z**



SCALER 2YS-3YS

AES2YS-3YS**X**

AES2YS-3YS**Z**



SCALER 204S

AES204S**X**

AES204S**Z**



SCALER 204S XP

XP

AES204SXP**X**

AES204SXP**Z**



SCALER 204SD

AES204SD**X**

AES204SD**Z**




SCALER 204SD XP

XP

AES204SDXP**X**

AES204SDXP**Z**



SCALER 311-312

AES311-312**X**



SCALER 311-312 XP

XP

AES311-312XP**X**

AES311-312XP**Z**



SCALER BATES 7-8

AESB7-8**X**

AESB7-8**Z**



SCALER BOGE 513*

XP

AESB513XP**X**

AESB513XP**Z**

**A end is a hoe, B end is a N128*



SCALER BLACKJACK*

AESBJ**X**

AESBJ**Z**

**4mm extended shank*



SCALER BLACKJACK XP*

XP THIN

AESBJXP**X**

AESBJXP**Z**

AESBJTXP**X**

**4mm extended shank*



SCALER EAGLE CLAW

AESEC**X**

AESEC**Z**



SCALER EAGLE CLAW XP

XP

AESECXP**X**


AESECXP**Z**



SCALER EAGLE TALON

AESET**X**

AESET**Z**




SCALER EAGLE TALON XP

XP

AESETXP**X**

AESETXP**Z**



SCALER H5-33

AESH5-33XP**X**

AESH5-33XP**Z**



SCALER H5-33 XP

XP

AESH5-33XPXP**X**

AESH5-33XPXP**Z**



SCALER H6-7

AESH6-7**X**

AESH6-7**Z**




SCALER H6-7 XP

XP

AESH6-7XP**X**

AESH6-7XP**Z**



SCALER J14-15

AESJ14-15**X**

AESJ14-15**Z**

AVAILABLE EXCLUSIVELY IN STAINLESS STEEL



SCALER J3-U15

AESJ3-U15**X**

AESJ3-U15**Z**



SCALER J3-U15 XP

XP

AESJ3-U15XP**X**

SCALERS *(continued)*



SCALER J3-U5

AESJ3-U5**X**




SCALER J3-U5 XP

XP

AESJ3-U5XP**X**

AESJ3-U5XP**Z**



SCALER J3-4

AESJ3-4**X**

AESJ3-4**Z**



SCALER J30-33

AESJ30-33**Z**



SCALER J31-32


AESJ31-32**Z**

AVAILABLE EXCLUSIVELY IN STAINLESS STEEL



SCALER J4-U5

AESJ4-U5**X**



SCALER J4-U5 XP

XP

AESJ4-U5XP**X**

AESJ4-U5XP**Z**



SCALER MORSE 0-00

AESM0-00**X**

AESM0-00**Z**



SCALER N135

AESN135**X**



SCALER M23

AESM23**X**

AESM23**Z**



SCALER M23XP

XP THIN

AESM23XP**X**

AESM23XP**Z**

AESM23TXP**X**

AESM23TXP**Z**



SCALER NEB 128 DOUBLE-END

AESN128DE**X**

AESN128DE**Z**



SCALER NEB 128 DOUBLE-END XP

XP

AESN128DEXP**X**

AESN128DEXP**Z**



SCALER N1

AESN1**X**

AESN1**Z**



SCALER N1 XP

XP

AESN1XP**X**

AESN1XP**Z**



SCALER N2

AESN2**X**

AESN2**Z**



SCALER N2 XP

XP

AESN2XP**X**

AESN2XP**Z**



SCALER N137

AESN137**X**



SCALER N137 XP

XP

AESN137XP**X**

AESN137XP**Z**



SCALER NEB 128 SINGLE-END

AESN128**X**



SCALER NEB 128 SINGLE-END XP

XP

AESN128XP**X**

AESN128XP**Z**




SCALER U15 SINGLE-END

AESU15**X**

AVAILABLE EXCLUSIVELY IN STAINLESS STEEL

SCALERS *(continued)*



SCALER U15-30


AESU15-30**X**



SCALER U15-30 XP

XP

AESU15-30XP**X**



SCALER U15-33

AESU15-33**X**

AESU15-33**Z**

AVAILABLE EXCLUSIVELY IN STAINLESS STEEL

KITS



Available only in resin handles.

KIT 1 - EAGLE CLAW™ & EAGLE TALON™ I AEGEKXPX

The universal scaler set— a combination of the Talon for anterior and the Claw for posterior.

Eagle Talon is a combination of two already popular anterior tips: the Jacquette 3 and the Nebraska 128.

Eagle Claw is the ideal universal posterior scaler with a classic design for tight interproximal areas.

XP® Sharpen-Free instruments keep their shape throughout the useful life of the instrument.



Available only in resin handles.

KIT 2 - GRACEY 11-12 & GRACEY 13-14 I AEGEK2XPX

The classic molar-mesial and molar-distal Graceys.

XP means no sharpening – these are the first Graceys that retain their original shape.

Slim design – ideal for heavy calculus and for tight interproximal areas of recall patients.



Available only in resin handles.

KIT 6 M23 & BLACKJACK I AEGEK6XPX

Blackjack - The ultimate anterior scaler.


M23 - Replaces sickles and mini-sickles for posterior.

Ideal for safe removal of deposits.

XP Sharpen-Free Technology means the instruments retain their original shape.

Extra-light handle design that helps reduce repetitive motion injuries.

IMPLANT INSTRUMENTS



IMPLANT INSTRUMENT 204S

AEII204S**X**



IMPLANT INSTRUMENT B5-6

AEIIB5-6**X**



IMPLANT INSTRUMENT L3-4

AEIIL3-4**Z**



IMPLANT INSTRUMENT N128-L5


AEIIN128-L5**X**

IMPLANT INSTRUMENT KITS



IMPLANT KIT WITH PLASTIC CASSETTE I AEIIKP

Both implant instrument kits include all four implant instruments seen above.



IMPLANT KIT WITH SL5 CASSETTE I AEIIKZ

Both implant instrument kits include all four implant instruments seen above.

FILES



HIRSCHFELD 3-7 FILE

AEFH3-7



HIRSCHFELD 5-11 FILE

AEFH5-11



1/4'' - 6.3MM **STAINLESS STEEL**

Both files compatible with 1/4 Stainless Steel handle.

QUIK-TIP™ CUSTOMIZATION

27 QUIK-TIP™ HANDLES AND WRENCH

28 DOUBLE GRACEY

28 STANDARD GRACEY

30 GRACEY ACCESS

31 GRACEY DEEP POCKET

32 CURETTES

33 SCALERS

QUIK-TIP™ CUSTOMIZATION

Designed with the hygienist in mind, Quik-Tips offer the ability to discard a single end when it becomes dull without having to replace your entire instrument, **giving your practice significant cost savings.**

Quik-Tip cone socket benefits:

- Completely customizable
- Significant cost savings
- Safe alternative to re-tipping
- Eight handle colors to choose from
- Re-tip in seconds
- Reduce waste



QUIK-TIP HANDLES



DOUBLE-ENDED HANDLE, BLUE I AEQTHB



DOUBLE-ENDED HANDLE, GREEN I AEQTHG



DOUBLE-ENDED HANDLE, ORANGE I AEQTHO



DOUBLE-ENDED HANDLE, PURPLE I AEQTHP



DOUBLE-ENDED HANDLE, RED I AEQTHR



DOUBLE-ENDED HANDLE, BLACK I AEQTHX



DOUBLE-ENDED HANDLE, YELLOW I AEQTHY



DOUBLE-ENDED HANDLE, STAINLESS STEEL I AEQTHZ



QUIK-TIP WRENCH I AEQTW

Remember to order a Quik-Tip wrench with your Quik-Tip order!

DOUBLE GRACEY™ - QUIK-TIP™

The Double Gracey line of instruments started with the Scandette (now the Double Gracey Posterior) and later transformed into a complete line of dual-bladed curettes. For more information on Double Gracey curettes, go to pages 16 and 17.

AVAILABLE EXCLUSIVELY IN XP SHARPEN-FREE



DOUBLE GRACEY ANTERIOR (A)

AEDGAAXPQT





DOUBLE GRACEY ANTERIOR (B)

AEDGABXPQT



DOUBLE GRACEY POSTERIOR (A)

AEDGPAXPQT



DOUBLE GRACEY POSTERIOR (B)

AEDGPBXPQT



MINI DOUBLE GRACEY ANTERIOR (A)

AEDGMAAXPQT



MINI DOUBLE GRACEY ANTERIOR (B)

AEDGMABXPQT



MINI DOUBLE GRACEY POSTERIOR (A)

AEDGMPAXPQT



MINI DOUBLE GRACEY POSTERIOR (B)

AEDGMPBXPQT

STANDARD GRACEY - QUIK-TIP



GRACEY 1 XP

AEG1XPQT XP



GRACEY 2 XP

AEG2XPQT XP



GRACEY 3 XP

AEG3XPQT XP



GRACEY 4 XP

AEG4XPQT XP



GRACEY 5 XP

AEG5XPQT XP



GRACEY 6 XP

AEG6XPQT XP



GRACEY 7 XP

AEG7XPQT XP



GRACEY 8 XP

AEG8XPQT XP

STANDARD GRACEY - QUIK-TIP™ *(continued)*



GRACEY 9 XP

XP

AEG9XPQT



GRACEY 10 XP

XP

AEG10XPQT



GRACEY 11 XP

XP

AEG11XPQT



GRACEY 12 XP

XP

AEG12XPQT



GRACEY 13 XP

XP

AEG13XPQT



GRACEY 14 XP

XP

AEG14XPQT



GRACEY 15 XP

XP

AEG15XPQT



GRACEY 16 XP

XP

AEG16XPQT



GRACEY 17 XP

XP

AEG17XPQT



GRACEY 18 XP

XP

AEG18XPQT



GRACEY 1

AEG1TTQT



GRACEY 2

AEG2TTQT



GRACEY 5

AEG5TTQT



GRACEY 6

AEG6TTQT



GRACEY 7

AEG7TTQT



GRACEY 8

AEG8TTQT



GRACEY 11

AEG11TTQT



GRACEY 12

AEG12TTQT



GRACEY 13

AEG13TTQT



GRACEY 14

AEG14TTQT

GRACEY ACCESS - QUIK-TIP™ *(50% Shorter Blade and Extended Shank)*



GRACEY ACCESS 0 XP

XP

AEGA00XPQT



GRACEY ACCESS 00 XP

XP

AEGA0XPQT



GRACEY ACCESS 1 XP

XP THIN

AEGA1XPQT

AEGA1TXPQT



GRACEY ACCESS 2 XP

XP THIN

AEGA2XPQT

AEGA2TXPQT



GRACEY ACCESS 3 XP

XP

AEGA3XPQT



GRACEY ACCESS 4 XP

XP

AEGA4XPQT



GRACEY ACCESS 5 XP

XP

AEGA5XPQT



GRACEY ACCESS 6 XP

XP

AEGA6XPQT



GRACEY ACCESS 7 XP

XP THIN

AEGA7XPQT

AEGA7TXPQT



GRACEY ACCESS 8 XP

XP THIN

AEGA8XPQT

AEGA8XPQT



GRACEY ACCESS 11 XP

XP

AEGA11XPQT



GRACEY ACCESS 12 XP

XP

AEGA12XPQT



GRACEY ACCESS 13 XP

XP THIN

AEGA13XPQT

AEGA13TXPQT



GRACEY ACCESS 14 XP

XP THIN

AEGA14XPQT

AEGA14TXPQT



GRACEY ACCESS 15 XP

XP THIN

AEGA15XPQT

AEGA15TXPQT



GRACEY ACCESS 16 XP

XP THIN

AEGA16XPQT

AEGA15TXPQT

GRACEY DEEP POCKET - QUIK-TIP™



GRACEY DEEP POCKET 1 XP

XP

AEGDP1XPQT



GRACEY DEEP POCKET 2 XP

XP

AEGDP2XPQT



GRACEY DEEP POCKET 4 XP

XP

AEGDP4XPQT



GRACEY DEEP POCKET 5 XP

XP

AEGDP5XPQT



GRACEY DEEP POCKET 7 XP

XP

AEGDP7XPQT



GRACEY DEEP POCKET 8 XP

XP

AEGDP8XPQT



GRACEY DEEP POCKET 11 XP

XP

AEGDP11XPQT



GRACEY DEEP POCKET 12 XP

XP

AEGDP12XPQT



GRACEY DEEP POCKET 13 XP

XP

AEGDP13XPQT



GRACEY DEEP POCKET 14 XP

XP

AEGDP14XPQT



GRACEY DEEP POCKET 15 XP

XP

AEGDP15XPQT



GRACEY DEEP POCKET 16 XP

XP

AEGDP16XPQT

CURETTES - QUIK-TIP™



BARNHART 5

AECB5TTQT



BARNHART 5 XP

XP THIN

AECB5XPQT

AECB5TXPQT



BARNHART 6

AECB6TTQT



BARNHART 6 XP

XP THIN

AECB6XPQT

AECB6TXPQT



COLUMBIA 13

AECC13TTQT



COLUMBIA 13 XP

XP

AECC13XPQT



COLUMBIA 14

AECC14TTQT



COLUMBIA 14 XP

XP

AECC14XPQT



COLUMBIA 4L

AECC4LTTQT



COLUMBIA 4L XP

XP

AECC4LXPQT



COLUMBIA 4R

AECC4RTTQT



COLUMBIA 4R XP

XP

AECC4RXPQT



LANGER 1 XP

XP

AECL1XPQT



LANGER 2 XP

XP

AECL2XPQT



LANGER 3 XP

XP

AECL3XPQT



LANGER 4 XP

XP

AECL4XPQT



LANGER 5 XP

XP

AECL5XPQT



LANGER 6 XP

XP

AECL6XPQT



CURETTE 411 XP

XP

AEC411XPQT



CURETTE 412 XP

XP

AEC412XPQT



YOUNGER-GOOD 7

AECY7TTQT



YOUNGER-GOOD 7 XP

XP

AECY7XPQT



YOUNGER-GOOD 8

AECY8TTQT



YOUNGER-GOOD 8 XP

XP

AECY8XPQT

1-800-551-5172 | Fax: 406-549-7452 | www.am-eagle.com

CURETTES - QUIK-TIP™ *(continued)*



RULES 3

AECR3TTQT



RULES 3 XP

XP

AECR3XPQT



RULES 4

AECR4TTQT



RULES 4 XP

XP

AECR4XPQT

SCALERS - QUIK-TIP



204S A

AES204SATTQT



204S A XP

XP

AES204SAXPQT



204S B

AES204SBTTQT



204S B XP

XP

AES204SBXPQT



204SD A

AES204SDATTQT



204SD A XP

XP

AES204SDAXPQT



204SD B

AES204SDBTTQT



204SD B XP

XP

AES204SDBXPQT



311

AES311TTQT



311 XP

XP

AES311XPQT



312

AES312TTQT



312 XP

XP

AES312XPQT



33

AES33TTQT



33 XP

XP

AES33XPQT



BOGE 513 A XP

XP

AESB513AXPQT



BOGE 513 B XP

XP

AESB513BXPQT



BLACKJACK A

AESBJATTQT



BLACKJACK A XP

XP THIN

AESBJAXPQT

AESBJTAXPQT



BLACKJACK B

AESBJBTTQT



BLACKJACK B XP

XP THIN

AESBJBXPQT

AESBJTBXPQT

SCALERS - QUIK-TIP™ *(continued)*



EAGLE CLAW A

AESECATTQT



EAGLE CLAW A XP

XP

AESECAXPQT



EAGLE CLAW B

AESECBTTQT



EAGLE CLAW B XP

XP

AESECBXPQT



EAGLE TALON A

AESETATTQT



EAGLE TALON A XP

XP

AESETAXPQT



EAGLE TALON B

AESETBTTQT



EAGLE TALON B XP

XP

AESETBXPQT



H5

AESH5TTQT



H5 XP

XP

AESH5XPQT



H6

AESH6TTQT



H6 XP

XP

AESH6XPQT



H7

AESH7TTQT



H7 XP

XP

AESH7XPQT



M23 A

AESM23ATTQT



M23 A XP

XP THIN

AESM23AXPQT

AESM23TAXPQT



M23 B

AESM23BTTQT



M23 B XP

XP THIN

AESM23BXPQT

AESM23TBXPQT



N1 A XP

XP

AESN1AXPQT



N1 B XP

XP

AESN1BXPQT



N128

AESN128TTQT



N128 XP

XP

AESN128XPQT

1-800-551-5172 | Fax: 406-549-7452 | www.am-eagle.com

DIAGNOSTIC

36 MIRRORS AND HANDLES

37 EXPLORERS

37 PROBES

38 EXPLORERS/PROBES

MIRRORS & HANDLES



MIRROR #4

AEM4

MIRROR #4 12 PACK

AEM4/12



MIRROR #5

AEM5

MIRROR #5 12 PACK

AEM5/12



1/4" - 6.3 MM ROUND I AEMHR *(shown)*



1/4" - 6.3 MM TAPERED I AEMHT *(shown)*



3/8" - 9.5 MM ROUND I AEMHRZ



3/8" - 9.5 MM BLACK ROUND I AEMHRX *(shown)*

3/8" - 9.5 MM BLACK SIMPLE STEM ROUND I AEMHRSSX



3/8" - 9.5 MM BLACK TAPERED I AEMHTX *(shown)*

3/8" - 9.5 MM BLACK SIMPLE STEM TAPERED I AEMHTSSX

EXPLORERS



EXPLORER DE 11-12

AEEXP11-12**X**
AEEXP11-12



EXPLORER 11-12 DEEP POCKET

AEEXP11-12DP**X**
AEEXP11-12DP



EXPLORER PIGTAIL 2

AEEXP2**X**
AEEXP2



EXPLORER SHEPHERD'S HOOK 23

AEEXP23**X**
AEEXP23



EXPLORER 23-6

AEEXP23-6**X**
AEEXP23-6



EXPLORER 3

XAEEXP3



EXPLORER 3 COW HORN

AEEXP3CH**X**
AEEXP3CH



EXPLORER 3ES

AEEXP3ES**X**
AEEXP3ES



EXPLORER 5

AEEXP5**X**
AEEXP5



EXPLORER 6 SE

AEEXP6SE**X**
AEEXP6SE

IMAGE NOT AVAILABLE

EXPLORER 9

AEEXP9**X**
AEEXP9

AVAILABLE INTERNATIONALLY ONLY

PROBES



PROBE MICHIGAN 0

AEP0B



PROBE MICHIGAN 0 WITH WILLIAMS MARKINGS

AEP0WB**X**
AEP0WB



PROBE WILLIAMS 1- GOLDMAN FOX 2

AEP1-2B**X**
AEP1-2B




PROBE 12

AEP12B**X**
AEP12B



PROBE 12

AEP12Y**X**
AEP12Y




PROBE 12

AEP12Y/G**X**
AEP12Y/G



PROBE WILLIAMS 1

AEP1B**X**
AEP1B



PROBE WILLIAMS 1

AEP1Y



PROBE 2-12

AEP2-12B**X**
AEP2-12B



PROBE 2-12

AEP2-12Y**X**
AEP2-12Y

PROBES *(continued)*



PROBE NABERS 2

AEPN2B

AEPN2B**X**



PROBE SOUTH DAKOTA 12

AEPSD12Y

AEPSD12Y**X**



PROBE UNC 15 RUNG

AEPUNC15RB

AEPUNC15RB**X**



PROBE BALL-TIP SCREENING

AEPWHOB


AEPWHOB**X**

EXPLORERS / PROBES



EXPLORER 23/ PROBE MICHIGAN 0


AEEP23/0B



EXPLORER 23/ PROBE 12

AEEP23/12B**X**

AEEP23/12B



EXPLORER 23/ PROBE 12

AEEP23/12Y**X**

AEEP23/12Y



EXPLORER 23/ PROBE 12

AEEP23/12YG**X**

AEEP23/12YG



EXPLORER 23/ PROBE WILLIAMS 1

AEEP23/1B**X**

AEEP23/1B



EXPLORER 23/ PROBE WILLIAMS 1

AEEP23/1Y**X**

AEEP23/1Y



EXPLORER 23/ SOUTH DAKOTA 12

AEEP23/SD12Y**X**

AEEP23/SD12Y



EXPLORER 23/ PROBE BALL-TIP SCREENING

AEEP23/WHOB**X**

AEEP23/WHOB

1-800-551-5172 | Fax: 406-549-7452 | www.am-eagle.com

RESTORATIVE

40 BURNISHERS

40 CARVERS

41 DISCOID-CLEOIDS

41 CONDENSERS / PLUGGERS

42 SUPER PLUGGERS

42 CHISELS

42 SPATULA

42 EXCAVATORS

43 HATCHETS

43 MARGIN TRIMMERS

44 PACKING INSTRUMENTS

44 COMPOSITE

45 DR. MOPPER COMPOSITE SERIES

45 DR. COLONNA COMPOSITE SERIES

46 DR. BELVEDERE COMPOSITE CONTACT SERIES

46 ENDO INSTRUMENTS

47 SPREADER INSTRUMENTS

47 WAXING INSTRUMENTS

BURNISHERS



BURNISHER PLACEMENT INSTRUMENT- SHOWN WITH STANDARD 1/4" - 6.3 MM STAINLESS STEEL DOUBLE END HANDLE

AEPPA



BURNISHER 2-29 BEAVERTAIL-FOOTBALL

AEB2-29BEF



BURNISHER 2-30 BEAVERTAIL-FOOTBALL

AEB2-30BEF



BURNISHER 21B ACORN

AEB21BA



BURNISHER 26-27 BALL

AEB26-27B



BURNISHER 26-28 BALL-FOOTBALL

AEB26-28BF



BURNISHER 26-29 BALL-FOOTBALL

AEB26-29BF



BURNISHER NYU B1

AEBNYUB-1



BURNISHER PFI 18 BALL

AEBPFI18B



BURNISHER PFI 21 ACORN

AEBPFI21A



BURNISHER PITTS PEDO 18

AEBPP18



LINING CEMENT/CALCIUM HYDROXIDE PLACEMENT INSTRUMENT

AEPBP-2

CARVERS



BURNISHER PLACEMENT INSTRUMENT- SHOWN WITH STANDARD 1/4" - 6.3 MM STAINLESS STEEL DOUBLE END HANDLE

AECAG4-5



CARVER PFI 13

AECA13



CARVER PFI 6

AECA6



CARVER PFI 8A

AECA8A



CARVER PFI D3

AECAD3



CARVER PFI ASH 6

AECAA-6



CARVER PFI S6

AECAS6

1-800-551-5172 | Fax: 406-549-7452 | www.am-eagle.com

CARVERS *(continued)*



CARVER HOLLENBACK 1/23

AECAH1/23



CARVER HOLLENBACK 3

AECAH3



CARVER TANNER 3

AECAT3



CARVER IPC

AECAIPC



CARVER MINNESOTA 1

AECAIPC1M



CARVER MINNESOTA 2

AECAIPC2M



CARVER GREGG 4-5 XP

AECAG4-5TN**Z**

DISCOID-CLEOIDS



DISCOID-CLEOID 1-2 - SHOWN WITH STANDARD 1/4" - 6.3 MM STAINLESS STEEL DOUBLE END HANDLE

AECADC1-2



DISCOID-CLEOID 3-6

AECADC3-6



DISCOID-CLEOID 4-5

AECADC4-5



DISCOID-CLEOID 5T

AECADC5T



DISCOID-CLEOID 89-92

AECADC89-92



DISCOID-CLEOID UW-D5

AECADCUW-D5

CONDENSERS / PLUGGERS



CONDENSER / PLUGGER 0 - SHOWN WITH STANDARD 1/4" - 6.3 MM STAINLESS STEEL DOUBLE END HANDLE

AEC/PO

CONDENSERS / PLUGGERS *(continued)*

1.4 mm — 2.0 mm

CONDENSER/ PLUGGER 1-2

AEC/P1-2

1.5 mm — 2.5 mm

CONDENSER/ PLUGGER 3-4

AEC/P3-4

1.0 mm — 1.5 mm

CONDENSER/ PLUGGER HOLLENBACK 1

AEC/PH1

0.8 mm — 1.0 mm

CONDENSER/ PLUGGER OREGON 1

AEC/POR1

1.0 mm — 1.5 mm

CONDENSER/ PLUGGER OREGON 2

AEC/POR2

1.3 mm — 2.3 mm

CONDENSER/ PLUGGER OREGON 3

AEC/POR3

SUPER PLUGGERS

SUPER PLUGGER SMALL/MEDIUM

AEC/PSPSM

SUPER PLUGGER MEDIUM/LARGE

AEC/PSPML

CHISELS

CHISEL WEDELSTEADT 5-6

AEC5-6W

CHISEL BINANGLE 11-12

AEC11-12B

SPATULA

CEMENT SPATULA 8

AECSP8

EXCAVATORS

1.2 mm

EXCAVATOR 17

AEEXC17

1.5 mm

EXCAVATOR 18

AEEXC18

1.15 mm

EXCAVATOR 38-39

AEEXC38-39

1-800-551-5172 | Fax: 406-549-7452 | www.am-eagle.com

EXCAVATORS *(continued)*



EXCAVATOR 40-41

AEEXC40-41



EXCAVATOR 42-43

AEEXC42-43



EXCAVATOR 64-65

AEEXC64-65



EXCAVATOR 66-67

AEEXC66-67



EXCAVATOR 70-71

AEEXC70-71



EXCAVATOR 74-75

AEEXC74-75



EXCAVATOR 78-79

AEEXC78-79



EXCAVATOR E-1

AEEXCE-1



EXCAVATOR E-2

AEEXCE-2



EXCAVATOR E-3

AEEXCE-3

HATCHETS



HATCHET 8-9

AEHA8-9



HATCHET 10-11

AEHA10-11

MARGIN TRIMMERS



MARGIN TRIMMER 26

AEMT26



MARGIN TRIMMER 27

AEMT27



MARGIN TRIMMER 28

AEMT28




MARGIN TRIMMER 29


AEMT29

DR. BLACK'S INSTRUMENT NUMBERING FORMULA

Dr. Black's formula is based on the concept of instrument balance. Balance is established by keeping the blade of the instrument within 2mm of the long axis of the handle line. Proper balance is necessary to achieve optimal cutting with minimal effort.



Some instruments, such as margin trimmers, are constructed with a cutting edge at an angle other than a right angle. In these cases, a fourth number is placed in the formula that indicates the degree of angle of the cutting edge to the long axis of the handle.



PACKING INSTRUMENTS



PACKING INSTRUMENT CSI-1 SERRATED

AEPCSI-1



PACKING INSTRUMENT GUYER 7 SERRATED

AEPG7



PACKING INSTRUMENT ROUND SERRATED

AEPRNDS

COMPOSITE



1/4" - 6.3 MM **STANDARD STAINLESS STEEL**



3/8" - 9.5MM EAGLELITE **STAINLESS STEEL (Z)**



COMPOSITE 1

AEC1TN**Z**



COMPOSITE 2

AEC2TN**Z**



COMPOSITE 3

AEC3TN**Z**



COMPOSITE 4

AEC4TN

AEC4TN**Z**



COMPOSITE 5

AEC5TN

AEC5TN**Z**



COMPOSITE 6-21 ACORN

AEC6-21ATN

AEC6-21ATN**Z**



COMPOSITE 6

AEC6TN

AEC6TN**Z**



COMPOSITE 7

AEC7TN

AEC7TN**Z**



COMPOSITE 9

AEC9TN

AEC9TN**Z**



CARVER GREGG 4-5

AECAG4-5TN

AECAG4-5TN**Z**

DR. MOPPER COMPOSITE SERIES



8A LONG CARVER

AECA8ALTN

AECA8ALTN**Z**



8A CARVER

AECA8ATN

AECA8ATN**Z**



INTERPROXIMAL LONG CARVER

AECAIPCLTN

AECAIPCLTN**Z**



INTERPROXIMAL THIN CARVER

AECAIPCTTN

AECAIPCTTN**Z**



SUPER PLUGGER M/L

AEC/PMSPMLTN


AEC/PMSPMLTN**Z**



SUPER PLUGGER S/M XP

AEC/PMSPSMTN


AEC/PMSPSMTN**Z**



26-30 BALL-FOOTBALL

AEB26-30BFTN

AEB26-30BFTN**Z**



M-1 BALL-BURNISHER

AEBM1TN

AEBM1TN**Z**

DR. COLONNA COMPOSITE SERIES



ACORN CARVER

AEMCB21ATN**Z**



CONDENSER/PLUGGER 0

AEMCC/P0TN**Z**



CONDENSER/PLUGGER 1

AEMCC/P1TN**Z**



INTERPROXIMAL CARVER

AEMCCAIPCTN**Z**



MICRO-SPOON EXCAVATOR 4

AEMCEXC4TN**Z**



MICRO-SPOON EXCAVATOR 6

AEMCEXC6TN**Z**

DR. BELVEDERE COMPOSITE CONTACT SERIES



5/16'' - 7.9 MM **STAINLESS STEEL BELVEDERE DOUBLE END**



CONTACT FORMER REGULAR

AECCF



CONTACT FORMER SMALL


AECCFSM



CARVER SMALL COMPACT FORMER


AECP1CFSM

ENDO INSTRUMENTS




ROOT CANAL 1-3

AEE1-3P




ROOT CANAL 3-5

AEE3-5P



ROOT CANAL 5-7

AEE5-7P



ROOT CANAL 9-11

AEE9-11P



ENDO EXCAVATOR 31L

AEEM31LEXC



ENDO EXPLORER DG 16

AEEMDG16



ENDO CARVER GLICK 1

AEEMG1CA

SPREADER INSTRUMENTS



SPREADER D11-21

AEED11-21S



SPREADER THIN D11T-21

AEED11T-21S



SPREADER THIN D11T-25

AEED11T-25S

WAXING INSTRUMENTS



P.K. THOMAS 1

AEWPKT1



P.K. THOMAS 2

AEWPKT2



P.K. THOMAS 3

AEWPKT3



WAX SPATULA 7

AEWS7

FILES



HIRSCHFELD 3-7 FILE

AEFH3-7

Buccal/Lingual



HIRSCHFELD 5-11 FILE

AEFH5-11

Mesial/Distal



1/4'' - 6.3MM **STAINLESS STEEL**

Both files compatible with 1/4 Stainless Steel handle.

AUXILIARY PRODUCTS

49 COTTON PLIERS

50 INSTRUMENT CASSETTES

51 SHARPENING STONES AND ACCESSORIES

COTTON PLIERS



FLAGG 17 PLIER I AEPL17



COLLEGE 16 PLIER I AEPL4078



COLLEGE 16 LOCKING PLIER I AEPL4088

INSTRUMENT CASSETTES



CASSETTE 108 I AEOVS108

Holds 8 instruments



CASSETTE 109/12P I AEOVS109/12P

Holds 12 instruments



CASSETTE 110 I AEOVS110

Holds 10 instruments



CASSETTE 114 I AEOVS114

Holds 14 instruments

INSTRUMENT CASSETTES *(continued)*



CASSETTE 5 I AEOVSSL5R

Holds 5 instruments



CASSETTE 7 I AEOVSSL7

Holds 7 instruments



CASSETTE 8 I AEOVSWV8/8


Holds 8 instruments PLUS access area

1-800-551-5172 | Fax: 406-549-7452 | www.am-eagle.com

SHARPENING STONES & ACCESSORIES



SIMPLY SHARP CERAMIC STONE I AESACC/F



Rounded Edge:
Use for Curettes

Notched Side:
Course Grit

Opposite Side:
Fine Grit

ALL YOUR SHARPENING NEEDS IN ONE EASY KIT

SIMPLY SHARP CERAMIC STONE KIT I AESACC/FK

Sharpening dental instruments can be inconvenient, tedious and challenging. However, it is important to maintain sharp edges on your curettes and scalers. Our Simply Sharp Kit is designed with your convenience in mind. The double-sided stone offers you the choice of an extra-fine surface on one side and coarse grit on the other. The kit also includes a powerful 10x loupe for examining your edges. A sharpening test stick is included to conveniently test the cutting edge.

SHARPENING TECHNIQUE

1. Using the 10x magnification loupe, inspect cutting edges for defects.
2. Pick your method, sharpen your instrument and test (using the test stick).
3. Inspect again by using the magnification loupe for proper angle, edge contour and sharpness (a sharp edge will not reflect light).
4. Repeat steps 2 and 3. Make adjustments as needed until the desired edge is achieved.
5. Always check your instruments to maintain sharpness. We recommend a new instrument be kept for comparison.

SHARPENING STONES & ACCESSORIES *(continued)*



CERAMIC SUPER FINE STONE I AESACCS



DIAMOND STONE I AESADS



SHARPENING TEST STICK - 6 PACK I AESASTP *(shown)* I **SHARPENING TEST STICK** I AESASTS

To maintain optimal cutting edges on your instruments, American Eagle Instruments recommends routine sharpness testing by use of a test stick. To test for sharpness, simply stroke the test stick as you would during use. If the instrument "grabs," it is sharp; if it does not, it should be sharpened or discarded.

DO NOT sharpen XP® Sharpen-Free instruments. Simply discard and replace.

NOTES

NOTES

American Eagle Instruments, Inc. stands behind our entire product offering 100%. If you do not believe that American Eagle Instruments are of the finest quality, the unsatisfactory product may be returned within 30 days for exchange or credit. **This is our unconditional guarantee.**



6575 Butler Creek Road • Missoula, MT 59808 • 800-551-5172 • Fax: 406-549-7452 • www.am-eagle.com